**SERVOX® | SERVONA**

Absaugkatheter

suction catheter

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Catéter de succión

usisni kateter

Sací katéter

DE

EN

FR

IT

ES

HR

SK

**Orbisana Healthcare GmbH**

Biberweg 24-26  
53842 Troisdorf | Germany  
Tel. +49 (0)2241 9322-0  
Fax +49 (0)2241 9322-271  
info@orbisana.de  
www.orbisana.de

CE 0197

MD

CE0197

DE	CE-Kennzeichen mit Kennnummer der Benannten Stelle	Hersteller	Nicht wiederverwenden	Herstellungsdatum
EN	CE mark with the number of the notified body	Manufacturer	Do not re-use	Date of manufacture
FR	Marquage CE avec numéro d'identification de l'organisme notifié	Fabricant	Ne pas réutiliser	Date de fabrication
IT	Marchio CE con codice di identificazione dell'ente notificato	Produttore	Non riutilizzare	Data di produzione
ES	Marca CE con el número de identificación del organismo notificado	Fabricante	No reutilizable	Fecha de fabricación
HR	Oznaka CE s identifikacijskim brojem imenovane službe	Proizvođač	Nemojte ponovno upotrebljavati	Datum proizvodnje
SK	Označenie CE s číslom oznamovacieho orgánu	Výrobca	Nepoužívajte znova	Dátum výroby

LOT

DE	Gebrauchsanweisung beachten	Verwendbar bis	Nicht verwenden, wenn die Verpackung beschädigt ist	Chargenbezeichnung
EN	Consult instructions for use	Use-by date	Do not use if package is damaged	Batch code
FR	Observer le mode d'emploi	Date de péremption	Ne pas utiliser si l'emballage est endommagé	Numéro de lot
IT	Attenersi alle istruzioni per l'uso	Data di scadenza	Non utilizzare se la confezione è danneggiata	Numero del lotto
ES	Observe las instrucciones de uso	Fecha de caducidad	No utilizar si el embalaje está dañado	Número de lote
HR	Pridržavajte se uputa za uporabu	Rok uporabe	Ne upotrebljavajte ako je pakiranje	Broj šarže
SK	Prosím, poraďte sa s návodom na použitie	Dátum spotreby	Nepoužívať, ak je obal poškodený	Kód šarže

STERILE

REF

PZN

DE	Sterilisiert mit Ethylenoxid	Artikelnummer	Nicht erneut sterilisieren	Pharmazentralnummer
EN	Sterilized using ethylene oxide	Catalogue number	Do not resterilize	Central pharmaceutical number
FR	Stérilisé à l'oxyde d'éthylène	Référence	Ne pas restériliser	Numéro central pharmaceutique
IT	Sterilizzato con ossido di etilene	Codice articolo	Non risterilizzare	Numero centrale della farmacia
ES	Producto esterilizado con óxido de etileno	Número de referencia	No reesterilice	Número central de farmacia
HR	Sterilizirano etilenoksidom	Broj artikla	Nemojte ponovno sterilizirati	Broj Pharma centrale
SK	Sterilizované oxidom etylénu	Katalógové číslo	Neznovusterilizovať	Centrálné farmaceutické číslo

**Symbolerklärung | Symbol explanation | Explication des symboles | Spiegazione dei simboli | Explicación de símbolos | Objašnjenje simbola | Vysvetlenie symbolov**

				
<b>DE</b>	Vor Sonnenlicht schützen	Sicherheitshinweise beachten	Trocken aufbewahren	Temperaturbegrenzung
<b>EN</b>	Keep away from sunlight	Note the safety instructions	Keep dry	Temperature limit
<b>FR</b>	Protéger du rayonnement solaire	Respecter les consignes de sécurité	Conserver au sec	Limitation de la température
<b>IT</b>	Proteggere dalla luce solare	Attenersi alle avvertenze di sicurezza	Conservare all'asciutto	Limitazione di temperatura
<b>ES</b>	Mantener alejado de la luz solar	Observe las instrucciones de seguridad	Mantener seco	Límite de temperatura
<b>HR</b>	Čuvajte od sunčevog svjetla	Pridržavajte se sigurnosnih uputa	Čuvajte na suhom mjestu	Ograničenje temperature
<b>SK</b>	Chránit pred slnečným žiarením	Dodržiavajte bezpečnostné pokyny	Uchovávať v suchu	Teplotný limit

**DE**

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## SERVOX® | SERVONA

### suction catheter

<b>Symbol explanation</b>	<b>2</b>
<b>1 Foreword</b>	<b>10</b>
1.1 Intended purpose	10
1.2 Indication	11
1.3 Contraindication	11
1.4 Side effects	11
<b>2 Safety instructions</b>	<b>11</b>
<b>3 Description</b>	<b>12</b>
<b>4 Use</b>	<b>13</b>
<b>5 Storage, transport, disposal</b>	<b>13</b>
<b>6 Disclaimer</b>	<b>13</b>

## 1 Foreword

This instruction manual is intended for patients, relatives and medical staff. The patient can use the SERVOX® | SERVONA suction catheter independently (without the support of relatives or carers), but must be given instructions on its use by their doctor or by medical staff beforehand. SERVOX® | SERVONA suction catheters can be used at home and in a clinical setting. They are suitable for all age groups according to the size selected by the patient's doctor.

SERVOX® | SERVONA suction catheters should be used in line with the descriptions in this instruction manual and in accordance with the recognised rules of medicine.

### Note

Read this instruction manual carefully. It contains important information on safe and proper use of the product.

SERVOX® | SERVONA suction catheters are assigned to batches for traceability purposes. The batch number of each suction catheter is indicated on the packaging label next to "LOT". This batch number allows each suction catheter to be identified.

If you have any questions or problems with the suction catheters, please contact your field representative or Orbisana Healthcare GmbH, referencing this batch number.

### 1.1 Intended purpose

The SERVOX® | SERVONA suction catheters are intended for suctioning secretions from the respiratory tract and trachea. They are suitable for nasal, oral and endotracheal suctioning.

## Note

Suction catheters are sterile disposable products and must not be resterilised or reused.

### 1.2 Indication

For use on ventilated or spontaneously breathing tracheotomised, laryngectomised or intubated patients or patients who are not tracheotomised, laryngectomised or intubated but require help with expectoration.

### 1.3 Contraindication

Do not use on patients with skull fractures or patients with known allergies to the material used.

Patients with an acute tendency to bleed due to e.g. anticoagulation treatment may only use the SERVOX® | SERVONA suction catheter after consulting their doctor.

### 1.4 Side effects

There are no known side effects provided the patient uses the product properly.

EN

## 2 Safety instructions

Safe and proper use of your SERVOX® | SERVONA suction catheter is important for your health. Ignoring the safety instructions can have significant effects on your health or even cause life-threatening injuries.

**In the following, warnings on the safe use of the product are indicated with the  symbol.**

**Danger:** The word "Danger" indicates a likely situation that, if not avoided, can result in grave injury or even death.

**Warning:** The word "Warning" indicates a situation that, if not avoided, can result in grave injury or even death.

**Please observe the following safety instructions:**

### Danger

**Risk of injury/ suffocation:** The tube adhering to the tracheal wall, careless suctioning, suctioning too deep or using damaged suction catheters can irritate or injure the mucosa. This can cause the mucosa to swell and the trachea to seal, which can lead to suffocation. The vagus nerve can be irritated. This can cause cardiac arrest.

- Never use a damaged suction catheter!
- Always check the suction depth during suctioning!
- Always check the suction strength during suctioning!
- Keep the suction catheter moving during suctioning!
- During suctioning, try not to touch the trachea!

### Warning

**Risk of infection:** The use of dirty suction catheters or unhygienic work practices can cause infections.

- Do not use the suction catheter after the expiry date!
- Do not use the suction catheter if the packaging is damaged!

- Do not use a used suction catheter!
- Always ensure that you follow hygienic work practices!
- Wear disposable gloves during suctioning!

**Danger of suffocation:** Selecting the incorrect diameter of the suction catheter can interrupt breathing or cause negative pressure in the lungs.

- Always check the size of the suction catheter!

#### Note

Please inform us or the responsible authorities if you experience a life-threatening illness or injury or permanent physical damage as a result of using our product, or require hospital treatment or have to extend hospital treatment due to using our product.



## 3 Description

SERVOX® | SERVONA suction catheters are free of latex and phthalate and are available in various lengths, sizes, hardness grades and models. All SERVOX® | SERVONA suction catheters have a rounded tip.


SERVOX® | SERVONA suction catheters are used together with a suction unit and are comprised of a flexible PVC tube with a connection piece to connect it to the suction unit. To check the insertion depth, a scale in 1-cm increments is marked on the tube.

The lateral eyes on the sides of the tube above the rounded distal tip of the suction catheter prevent it from adhering to the mucosa. The atraumatic tip also has a bulbous head on the distal end of the catheter in order to minimise contact with the mucosa. The various tip shapes differ depending on how they are inserted into the trachea.

#### Connections:

	
<p><b>Funnel:</b> The funnel connection must be connected to the suction tube of a suction unit with an additional fingertip control.</p>	<p><b>Fingertip:</b> The fingertip connection can be used to connect the suction catheter to the suction tube on the suction unit.</p>

#### Tips:

		
<p><b>Atraumatic tip:</b> Bulbous head, rounded edge, 4 small offset eyes</p>	<p><b>Standard tip:</b> Rounded edge, 2 large offset eyes or 2 small eyes opposite each other</p>	
<p>→ Insert <b>with</b> suction</p>	<p>→ Insert <b>without</b> suction</p>	

**Size indicated by the coloured marking on the connection pieces:**

Size in Charrière (CH)	CH 6	CH 8	CH 10	CH 12	CH 14	CH 16	CH 18
Outer diameter in mm	2.0	2.67	3.33	4.0	4.67	5.33	6.0
Colour	Light green	Light blue	Black	White	Green	Orange	Red

## 4 Use

**Note**

- Always insert the suction catheter with the standard tip without suction.
- Always insert the suction catheter with the atraumatic tip with suction.

- Disinfect hands and put on gloves
- Check that the size, length and model of the suction catheter is correct
- Check the sterile packaging for the expiry date and damage. Do not use if the suction catheter has expired or the packaging is damaged
- Open the sterile packaging, take out the connection piece and slide it on the adapter or the fingertip control of the suction tube
- Switch on the suction unit (read the suction unit instruction manual)
- Then remove the rest of the suction catheter from the packaging
  - ➔ Insert the suction catheter with the **standard tip** carefully **without suction!**
  - ➔ Insert the suction catheter with the **atraumatic tip** carefully **with suction!**
- Then slowly pull the catheter out under suction with a slight twisting motion
- Do not suction for more than 15 s
- Dispose of the suction catheter after use

EN

## 5 Storage, transport, disposal

Store or transport the SERVVOX® | SERVONA suction catheters in a clean, dry environment that is free of dust and protected from light at temperatures between -5°C and +35°C.

Protect from direct sunlight and heat.

The SERVVOX® | SERVONA suction catheters must be disposed of in line with national regulations for waste products. Suction catheters can generally be disposed of in household waste. The packaging can be recycled.

## 6 Disclaimer

In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product.

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