SERVOX® Tracheo Pad

Saugfähige Tracheokomresse
Absorbent Tracheal Compress
Cuscinetto tracheostomico

CE 0123 MD

T&L Co., Ltd.
70-17, Wonam-ro, Won-gok-myeon
Anseong-si, Gyeonggi-do
17554 Korea
Tel. +82 31 651 5220
Fax. +82 31 651 6212

Orbisana Healthcare GmbH
Biberweg 24-26
53842 Troisdorf Germany
Tel. +49 (0)2241 9322-0
Fax +49 (0)2241 9322-277
info@orbisana.de
www.orbisana.de

EC REP

KTR Europe GmbH
Mergenthalerallee 77,
65760 Eschborn, Germany

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SERVOX[®] Tracheo Pad

Saugfähige Tracheokompresse

1 Vorwort

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Absorbent Tracheal Compress

1 Preface

These instructions for use are intended for patients, relatives and medical professionals in the home and clinical setting.

Note

Read these instructions for use carefully. They provide a lot of important information for safe and proper use. Keep these instructions for future reference.

For traceability, SERVOX[®] Tracheo Pads are assigned a batch number. The number is indicated on the packaging following the 'LOT' marking.

1.1 Intended use

The SERVOX[®] Tracheo Pad absorbs tracheal secretions to reduce irritation caused by heavy secretions at the tracheostoma and cushions the cannula shield to prevent pressure points.

1.2 Indication

Patients using a tracheostomy tube.

1.3 Contraindication

Allergies against polyurethane.

1.4 Side effects

If the indication, contraindication, safety and application instructions are observed, no side effects are reported.

2 Safety instructions

Not following these safety instructions can result in serious health effects or even life-threatening situations. The manufacturer accepts no liability in the event of damage caused by not following the instructions or by manipulation of the product.

Risk of infection

- Do not reuse!
- Do not use if sterile packaging is damaged!
- Do not use after expiry date!
- Replace after 24 hours at the latest or more frequently, if necessary!
- In case of heavy secretion, the compress may swell, or the front side may become discoloured. If that happens, you must replace the compress!

EN

Note

We ask you to inform us and, if necessary, the competent authority if you suffer a life-threatening illness or injury or permanent physical damage in connection with the use of our product or if you have to be treated as an in-patient or prolong an in-patient treatment due to the use of our product.

3 Description

The SERVOX® Seal Pad is a sterile polyurethane foam tracheostomy compress. The underside is a light-coloured, non-adhesive layer that absorbs fluids from the tracheostoma. The SERVOX® Tracheo Pad can be slid over the tracheostomy tube through the slit, even when the cannula is in situ.

4 Application

- Clean the tracheostoma as usual and dry the surrounding skin carefully.
- Cut the SERVOX® Tracheo Pad to size if necessary.
- Slide the SERVOX® Tracheo Pad with the lighter, fluid-absorbent side facing the patient over the tracheostomy tube under the cannula shield.
- For removing, lift the cannula shield slightly and pull out the SERVOX® Tracheo Pad downwards and discard it

5 Storage, transport and disposal

Store or transport in a dry and dust-free environment at temperatures between 1 °C and 30 °C. Protect from direct sunlight and heat. If no disease-related instructions (e.g. MRSA, ORSA) have to be observed, dispose of the pad with household waste

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